




**GASTROENTEROLOGY**  
**DISPOSABLE ENDOSCOPIC ACCESSORIES**

[www.promostmedikal.com](http://www.promostmedikal.com)

## CONTENT



ENDOSKOPIC HEMOCLIP .....01




SCLEROTHERAPY INJECTION NEEDLE..... 03




POLIPECTOMY SNARE ..... 05



DILATION BALLOONS..... 07



STONE EXTRACTION BASKET. .... 09



STONE EXTRACTION BALLOON..... 11

HYDROPHILIC GUIDEWIRE). .... 13

CYTOLOGY BRUSH..... 15

BOUGIE DILATOR SET. .... 16

BIOPSY FORCEPS. .... 17

BIOPSY FORCEPS ALLIGATOR TEETH.....19


SNARE WITH NET..... 21

GRASPING FORCEPS. .... 23




## ENDOSKOPIC HEMOCLIP

Usage : For Endoscopic Hemostasis and clamping the wound in GI tract




One-one rotation, accurate positioning




## ENDOSCOPIK HEMOCLIP



One-one rotation, accurate positioning.



Close the clip, grab to tissue, ready to release.




Open clip fully, ready to grab the tissue

## Specifications

| Model             | Clip opening size (mm) | Working Length (mm) | Endoscope Channel (mm) | Characteristics<br>Karakter özel |         |
|-------------------|------------------------|---------------------|------------------------|----------------------------------|---------|
| PRM-CD-09-160-O   | 09                     | 1600                | ≥ 2.8                  | Gastroscopy                      | Uncoaed |
| PRM-CD-11-160-O ● | 11                     | 1600                | ≥ 2.8                  |                                  |         |
| PRM-CD-13-160-O   | 13                     | 1600                | ≥ 2.8                  |                                  |         |
| PRM-CD-09-230-O   | 09                     | 2300                | ≥ 2.8                  | Colonoscopy                      | Coating |
| PRM-CD-11-230-O ● | 11                     | 2300                | ≥ 2.8                  |                                  |         |
| PRM-CD-13-230-O   | 13                     | 2300                | ≥ 2.8                  |                                  |         |
| PRM-CD-09-160-P   | 09                     | 1600                | ≥ 2.8                  | Gastroscopy                      | Coating |
| PRM-CD-11-160-P   | 11                     | 1600                | ≥ 2.8                  |                                  |         |
| PRM-CD-13-160-P   | 13                     | 1600                | ≥ 2.8                  |                                  |         |
| PRM-CD-09-230-P   | 09                     | 2300                | ≥ 2.8                  | Colonoscopy                      | Coating |
| PRM-CD-11-230-P   | 11                     | 2300                | ≥ 2.8                  |                                  |         |
| PRM-CD-13-230-P   | 13                     | 2300                | ≥ 2.8                  |                                  |         |

● Common specifications

## SCLEROTHERAPY INJECTION NEEDLE




Endoscopic Injection treatment of esophageal and gastric varices. Endoscopic injection of submucosa in GI tract




Shield needles designed to prevent damage to working channel , Classic handle design and easy to control needle moving.


## SCLEROTHERAPY INJECTION NEEDLE



BEFORE INJECTION



INJECTION




AFTER INJECTION

### Specifications


| Model                  | O.D. (mm) | Working Length (mm) | Needle Diameter (G) (mm) | Needle Length (mm) | Endoscope Channel (mm) | Remark          |
|------------------------|-----------|---------------------|--------------------------|--------------------|------------------------|-----------------|
| PRM-IN-18-120-2304-A   |           | 1200                | 23(0.6)                  | 4                  | ≥ 2.0                  | Gastroscopy     |
| PRM-IN-18-120-2504-A   |           | 1200                | 25(0.5)                  | 4                  | ≥ 2.0                  |                 |
| PRM-IN-18-160-2304-A   |           | 1600                | 23(0.6)                  | 4                  | ≥ 2.0                  | Micro endoscopy |
| PRM-IN-18-160-2504-A   |           | 1600                | 25(0.5)                  | 4                  | ≥ 2.0                  |                 |
| PRM-IN-23-160-2104-A   |           | 1600                | 21(0.8)                  | 4                  | ≥ 2.8                  |                 |
| PRM-IN-23-160-2304-A   |           | 1600                | 23(0.6)                  | 4                  | ≥ 2.8                  |                 |
| PRM-IN-23-160-2504-A   |           | 1600                | 25(0.5)                  | 4                  | ≥ 2.8                  |                 |
| PRM-IN-23-160-2105-A   |           | 1600                | 21(0.8)                  | 5                  | ≥ 2.8                  |                 |
| ● PRM-IN-23-160-2305-A |           | 1600                | 23(0.6)                  | 5                  | ≥ 2.8                  |                 |
| ● PRM-IN-23-160-2505-A |           | 1600                | 25(0.5)                  | 5                  | ≥ 2.8                  |                 |
| PRM-IN-23-230-2304-A   |           | 2300                | 23(0.6)                  | 4                  | ≥ 2.8                  |                 |
| PRM-IN-23-230-2504-A   |           | 2300                | 25(0.5)                  | 4                  | ≥ 2.8                  |                 |

● Common specifications

## POLYPECTOMY SNARE




**Usage:** For excising polyps and other redundant tissues in GI tract, by high frequency electricity in combination with endoscope.




Superior one-one rotation, accurate positioning

## POLYPECTOMY SNARE


Easier to fit lesion, perfect fit without curl .  
Easy for cutting, excellent performance at blood coagulation, minimal smoke and crust Made of imported high-strength medical stainless steel, with good electrification.



Hexagonal



Crescent



Oval Model


### Specifications

| Model                  | Loop Width (mm) | Length (mm) | Endoscope Channel (mm) | Loop Shape (mm) | Characteristics            |
|------------------------|-----------------|-------------|------------------------|-----------------|----------------------------|
| PRM -SD-18-120-20-A1   | 20              | 1200        | ≥ 2.0                  | Oval            | Gastrointestinal endoscopy |
| PRM -SD-23-160-25-A1   | 25              | 1600        | ≥ 2.8                  |                 |                            |
| ● PRM -SD-23-160-30-A1 | 30              | 1600        | ≥ 2.8                  |                 |                            |
| PRM -SD-23-160-35-A1   | 35              | 1600        | ≥ 2.8                  |                 |                            |
| ● PRM -SD-23-230-30-A1 | 30              | 2300        | ≥ 2.8                  |                 |                            |
| PRM -SD-23-230-35-A1   | 35              | 2300        | ≥ 2.8                  |                 |                            |
| PRM -SD-23-160-15-B1   | 15              | 1600        | ≥ 2.8                  | Crescent        | Gastrointestinal endoscopy |
| PRM -SD-23-160-20-B1   | 20              | 1600        | ≥ 2.8                  | Hexagonal       | Gastrointestinal endoscopy |
| PRM -SD-23-160-25-C1   | 25              | 1600        | ≥ 2.8                  |                 |                            |
| PRM -SD-23-160-30-C1   | 30              | 1600        | ≥ 2.8                  |                 |                            |




● Common specifications



## DILATION BALLOONS



Usage :For the dilation for digestive and respiratory tract stenosis.



High pressure-resistant materials ensure safe dilation. Cylindrical balloon is not easy to deform, and after dilation, it is convenient to observe. Various balloon size for choose, meets different clinical demand

## DILATION BALLOONS

### Dilation balloon for ERCP

| Model            | Balloon O.D. (mm) | Balloon | Catheter Length (mm) | Working Pressure (atm) | Remak                                                                             |
|------------------|-------------------|---------|----------------------|------------------------|-----------------------------------------------------------------------------------|
| PRM-BD-06-40-200 | 6                 |         | 2000                 | 5                      | Suitable for Endoscopic channel $\geq 2.8\text{mm}$ ; Guide wire diameter: 0.035" |
| PRM-BD-08-40-200 | 8                 |         | 2000                 | 5                      |                                                                                   |
| PRM-BD-10-40-200 | 10                |         | 2000                 | 5                      |                                                                                   |
| PRM-BD-12-40-200 | 12                |         | 2000                 | 5                      |                                                                                   |
| PRM-BD-14-40-200 | 14                |         | 2000                 | 5                      |                                                                                   |
| PRM-BD-16-40-200 | 16                |         | 2000                 | 5                      |                                                                                   |

### Dilation balloon for PTCS

| Model           | Balloon O.D. (mm) | Balloon | Catheter Length (mm) | Working Pressure (atm) | Remak                                                                             |
|-----------------|-------------------|---------|----------------------|------------------------|-----------------------------------------------------------------------------------|
| PRM-BD-06-40-90 | 6                 |         | 900                  | 5                      | Suitable for Endoscopic channel $\geq 2.0\text{mm}$ ; Guide wire diameter: 0.035" |
| PRM-BD-08-40-90 | 8                 |         | 900                  | 5                      |                                                                                   |
| PRM-BD-10-40-90 | 10                |         | 900                  | 5                      | Endoscopic channel $\geq 2.8\text{mm}$                                            |

### Dilation balloon for Esophagus, Pylorus and Intestinal Tract

| Model            | Balloon O.D. (mm) | Balloon | Catheter Length (mm) | Working Pressure (atm) | Remak                                                                             |
|------------------|-------------------|---------|----------------------|------------------------|-----------------------------------------------------------------------------------|
| PRM-BD-10-40-180 | 10                |         | 1800                 | 5                      | Suitable for Endoscopic channel $\geq 2.0\text{mm}$ ; Guide wire diameter: 0.035" |
| PRM-BD-12-80-180 | 12                |         | 1800                 | 5                      |                                                                                   |
| PRM-BD-14-80-180 | 14                |         | 1800                 | 5                      |                                                                                   |
| PRM-BD-16-60-180 | 16                |         | 1800                 | 5                      |                                                                                   |
| PRM-BD-18-60-180 | 18                |         | 1800                 | 4                      |                                                                                   |
| PRM-BD-18-80-180 | 18                |         | 1800                 | 4                      | Endoscopic channel $\geq 3.2\text{mm}$ ; Not suitable for Guide wire.             |
| PRM-BD-20-60-180 | 20                |         | 1800                 | 4                      |                                                                                   |

### Dilation balloon for Cardia

| Model           | Balloon O.D. (mm) | Balloon Length (mm) | Catheter Length (mm) | Working Pressure (atm) | Remak                                            |
|-----------------|-------------------|---------------------|----------------------|------------------------|--------------------------------------------------|
| PRM-BD-25-80-90 | 25                | 80                  | 900                  | 2                      | Not suitable for Endoscope use                   |
| PRM-BD-30-80-90 | 30                | 80                  | 900                  | 2                      |                                                  |
| PRM-BD-35-80-90 | 35                |                     | 900                  | 2                      | Guidewire with diameter 0.035" can pass through. |



### Other sizes

| Model               | Balloon O.D. (mm) | Balloon Length (mm) | Catheter Length (mm) | Working Pressure (atm) | Remak                                                                             |
|---------------------|-------------------|---------------------|----------------------|------------------------|-----------------------------------------------------------------------------------|
| PRM-BD-06-30-110/90 | 6                 | 30                  | 1100/900             | 5                      | Suitable for Endoscopic channel $\geq 2.0\text{mm}$ ; Guide wire diameter: 0.032" |
| PRM-BD-08-30-110/90 | 8                 | 30                  | 1100/900             | 5                      |                                                                                   |
| PRM-BD-10-40-110/90 | 10                | 40                  | 1100/900             | 5                      | Suitable for Endoscopic channel $\geq 2.8\text{mm}$ ; Guide wire diameter: 0.035" |
| PRM-BD-12-40-110/90 | 12                | 40                  | 1100/900             | 5                      |                                                                                   |
| PRM-BD-14-60-110/90 | 14                | 60                  | 1100/900             | 5                      |                                                                                   |
| PRM-BD-16-60-110/90 | 16                | 60                  | 1100/900             | 5                      |                                                                                   |
| PRM-BD-18-60-110/90 | 18                | 60                  | 1100/900             | 4                      |                                                                                   |




## STONE EXTRACTION BASKET

Remove gallstone in biliary duct or foreign bodies in upper and lower digestive tract.




## STONE EXTRACTION BASKET


Innovational handle design, with functions of push, pull and rotation, easier to grasp gallstone and foreign body. Convenient to inject contrast medium with injection port on handle.  
Made by advanced alloyed materials, ensures good shape retention even after difficult stone removal. Accept customization, can meet different requirements.




Diamond shape



Spiral shape



Tipless shape




Semi-circle shape


### Specifications

| Model                  | O.D. (mm) | Lenght(mm) | Basket Width (mm) | Endoscope Channel (mm) | Characteristics   |
|------------------------|-----------|------------|-------------------|------------------------|-------------------|
| PRM-BAS-18-70-15-S4-A  |           |            | 15                | ≥ 2.0                  | Diamond shape     |
| PRM-BAS-18-70-20-S4-A  |           |            | 20                | ≥ 2.0                  |                   |
| PRM-BAS-23-200-25-S4-A |           | 2000       | 25                | ≥ 2.8                  |                   |
| PRM-BAS-23-200-30-S4-A |           | 2000       | 30                | ≥ 2.8                  |                   |
| PRM-BAS-23-200-35-S4-A |           | 2000       | 35                | ≥ 2.8                  |                   |
| PRM-BAS-18-70-15-N4-B  |           |            | 15                | ≥ 2.0                  | Spiral shape      |
| PRM-BAS-18-70-20-N4-B  |           |            | 20                | ≥ 2.0                  |                   |
| PRM-BAS-18-70-15-N4-C  |           |            | 15                | ≥ 2.0                  |                   |
| PRM-BAS-18-70-20-N4-C  |           |            | 20                | ≥ 2.0                  | Semi-circle shape |
| PRM-BAS-23-200-20-N4-C |           | 2000       | 20                | ≥ 2.8                  |                   |
| PRM-BAS-23-200-25-N4-C |           | 2000       | 25                | ≥ 2.8                  |                   |
| PRM-BAS-23-200-30-N4-C |           | 2000       | 30                | ≥ 2.8                  |                   |
| PRM-BAS-18-70-15-N4-D  |           |            | 15                | ≥ 2.0                  | Tipless shape     |
| PRM-BAS-18-70-20-N4-D  |           |            | 20                | ≥ 2.0                  |                   |
| PRM-BAS-23-200-20-N4-D |           | 2000       | 20                | ≥ 2.8                  |                   |
| PRM-BAS-23-200-25-N4-D |           | 2000       | 25                | ≥ 2.8                  |                   |
| PRM-BAS-23-200-30-N4-D |           | 2000       | 30                | ≥ 2.8                  |                   |

## STONE EXTRACTION BALLOON




Stone removal in biliary tract, including biliary sludge and residual gallstone left in biliary tract after Mechanical lithotripsy




## STONE EXTRACTION BALLOON



Innovational design of three lumens. (lumens for Guide wire, Injection and Inflation), makes the insertion and radiography easier and more convenient



Unique design of tapering catheter improves the insertion into papilla. Imported balloon with a blast resistance which can help improve the safety of dilation.



Equipped with radiopaque marks on both ends of balloon, convenient for observing under X-Ray. Expand the channel for contrast medium going through more smoothly

### Specifications



| Model                    | Endoscope Channel (mm) | O.D. (mm)   | Length (mm) | Balloon Size (mm) |
|--------------------------|------------------------|-------------|-------------|-------------------|
| PRM-BAL-23-8.5/11.5/15-B | ≥ 3.2                  | 8.5/11.5/15 | 2000        | 7.5Fr             |
| PRM -BAL-23-9/12/15-B    | ≥ 3.2                  | 9/12/15     | 2000        | 7.5Fr             |
| PRM-BAL-23-10/13/16-B    | ≥ 3.2                  | 10/13/16    | 2000        | 7.5Fr             |

Common specifications for ERCP






## HYDROPHILIC GUIDE WIRE

To guide insertion of catheters and etc in non-vascular.



## HYDROPHILIC GUIDE WIRE



Anti-twist inner Niti core wire, offering an excellent twisting and pushing force. Super smooth PTFE zebra coating, easier to pass through the working channel, without any stimulation for tissue. Yellow & black coating, easier to track the guide wire and obvious under X-ray. Innovative triple anti-drop design at Hydrophilic tip, with no risk of drop-off. Straight tip design and angled tip design, providing more control options for doctors. We can accept customized service, such as the blue and white coating.


### Specifications

| Model            | O.D. (mm)      | Lenght(mm) | Characteristics                    |
|------------------|----------------|------------|------------------------------------|
| PRM-GW-80-150-C1 | 0.80 (0.032'') |            | C1: straight hydrophilic tip, 5 cm |
| PRM-GW-80-260-C1 | 0.80 (0.032'') |            |                                    |
| PRM-GW-80-450-C1 | 0.80 (0.032'') |            |                                    |
| PRM-GW-88-150-C1 | 0.88 (0.035'') |            |                                    |
| PRM-GW-88-260-C1 | 0.88 (0.035'') |            |                                    |
| PRM-GW-88-450-C1 | 0.88 (0.035'') |            |                                    |
| PRM-GW-63-150-C1 | 0.63 (0.025'') |            |                                    |
| PRM-GW-63-260-C1 | 0.63 (0.025'') |            |                                    |
| PRM-GW-63-450-C1 | 0.63 (0.025'') |            |                                    |
| PRM-GW-63-150-C4 | 0.63 (0.025'') |            | C4: angled hydrophilic tip, 5 cm   |
| PRM-GW-63-450-C4 | 0.63 (0.025'') |            |                                    |
| PRM-GW-80-450-C4 | 0.80 (0.032'') |            |                                    |
| PRM-GW-88-150-C4 | 0.88 (0.035'') |            |                                    |
| PRM-GW-88-450-C4 | 0.88 (0.035'') |            |                                    |
| PRM-GW-63-150-E1 | 0.63 (0.025'') |            | E1: straight hydrophilic tip, 9 cm |
| PRM-GW-63-450-E1 | 0.63 (0.025'') |            |                                    |
| PRM-GW-80-450-E1 | 0.80 (0.032'') |            |                                    |
| PRM-GW-88-150-E1 | 0.88 (0.035'') |            |                                    |
| PRM-GW-88-450-E1 | 0.88 (0.035'') |            |                                    |
| PRM-GW-63-150-E4 | 0.63 (0.025'') |            | E4: angled hydrophilic tip, 9 cm   |
| PRM-GW-63-450-E4 | 0.63 (0.025'') |            |                                    |
| PRM-GW-80-450-E4 | 0.80 (0.032'') |            |                                    |
| PRM-GW-88-150-E4 | 0.88 (0.035'') |            |                                    |
| PRM-GW-88-450-E4 | 0.88 (0.035'') |            |                                    |



CYTOLGY BRUSH

It is intended to be used for sampling tissues in clinical application.



Innovational brush design, with no risk of drop-off. Straight-shaped brush: easy to enter the depths of respiratory and digestive tract. U-Shaped brush: obtain tissues from all direction with a large quantity. Excellent price-performance-ratio.

Specifications

| Model               | CATHETER<br>O.D | Endoscope<br>Channel (mm) | Working<br>Length (mm) | Brush<br>O.D (mm) | Brush<br>Length( mm) | Characteristics             |
|---------------------|-----------------|---------------------------|------------------------|-------------------|----------------------|-----------------------------|
| PRM-CY-BC-18-120-A3 | 1.8             | 2.0                       | 1200                   | 3                 | 10                   | A: Straight Shaped<br>Brush |
| PRM-CY-BC-23-160-A3 | 2.3             | 2.8                       | 1600                   | 3                 | 8                    |                             |
| PRM-CY-BC-23-230-A3 | 2.3             | 2.8                       | 2300                   | 3                 | 8                    |                             |
| PRM-CY-BC-18-120-B  | 1.8             | 2.0                       | 1200                   | 2                 | 10                   | B: U Type Shape<br>Brush    |
| PRM-CY-BC-23-160-B  | 2.3             | 2.8                       | 1600                   | 3                 | 8                    |                             |
| PRM-CY-BC-23-230-B  | 2.3             | 2.8                       | 2300                   | 3                 | 8                    |                             |

## Buji Dilatör Set

**Intended use** Used for dilatation of esophageal stenosis.


Made by medical macromolecule, with proper rigidity.

Equipped with 2 radiopaque marks on each bougie, easy to observe under X-ray.

Various bougie sizes, meets different requirements.

Each set equipped with 6 bougies, 5, 7, 9, 13, 15, mm) 1 guidewire, 1 reusable cleaning brush and 1 flusher syringe.

With graduations at both sides, improved smooth tapered tip, no harm to guidewires or tissues.




### Specifications

| Model | Balloon Size (mm) | Length (mm) | Characteristics                                                                                 |
|-------|-------------------|-------------|-------------------------------------------------------------------------------------------------|
| PRM   | 5,7,9,11,13,15    | 850         | Each set equipped with 6 bougies, 1 guidewire, 1 reusable cleaning brush and 1 flusher syringe. |

## DISPOSABLE BIOPSY FORCEPS

Tissue sampling in digestive and respiratory tracts



## DISPOSABLE BIOPSY FORCEPS

Special sharpness treatment on cups, easier to sample. Innovative 4-bar linkage design, offering an excellent rigidity. Tiny interval between jaws, large occlusal force and outstanding fitness. Smooth jaws, effectively decrease the friction to working channel. Standard cups design with proper sampling capacity, small tissue harm.

### Disposable Biopsy Forceps (without coating) Specifications

| Model             | O.D. | Length (mm) | Endoscope Channel (mm) | Characteristics         |
|-------------------|------|-------------|------------------------|-------------------------|
| PRM-FB-18-120-O-O | 1.8  |             | ≥2.0                   | without spike, uncoated |
| PRM-FB-18-120-W-O | 1.8  |             | ≥2.0                   | with spike, uncoated    |
| PRM-FB-18-160-O-O | 1.8  |             | ≥2.0                   | without spike, uncoated |
| PRM-FB-18-160-W-O | 1.8  |             | ≥2.0                   | with spike, uncoated    |
| PRM-FB-23-160-O-O | 2.3  |             | ≥2.8                   | without spike, uncoated |
| PRM-FB-23-160-W-O | 2.3  |             | ≥2.8                   | with spike, uncoated    |
| PRM-FB-23-180-O-O | 2.3  |             | ≥2.8                   | without spike, uncoated |
| PRM-FB-23-180-W-O | 2.3  |             | ≥2.8                   | with spike, uncoated    |
| PRM-FB-23-200-O-O | 2.3  |             | ≥2.8                   | without spike, uncoated |
| PRM-FB-23-200-W-O | 2.3  |             | ≥2.8                   | with spike, uncoated    |
| PRM-FB-23-230-O-O | 2.3  |             | ≥2.8                   | without spike, uncoated |
| PRM-FB-23-230-W-O | 2.3  |             | ≥2.8                   | with spike, uncoated    |

### Disposable Biopsy Forceps (with coating) Specifications

| Model             | O.D. |  | Endoscope Channel (mm) | Characteristics       |
|-------------------|------|--|------------------------|-----------------------|
| PRM-FB-18-120-O-P | 1.8  |  | ≥2.0                   | without spike, coated |
| PRM-FB-18-120-W-P | 1.8  |  | ≥2.0                   | with spike, coated    |
| PRM-FB-18-160-O-P | 1.8  |  | ≥2.0                   | without spike, coated |
| PRM-FB-18-160-W-P | 1.8  |  | ≥2.0                   | with spike, coated    |
| PRM-FB-23-160-O-P | 2.3  |  | ≥2.8                   | without spike, coated |
| PRM-FB-23-160-W-P | 2.3  |  | ≥2.8                   | with spike, coated    |
| PRM-FB-23-180-O-P | 2.3  |  | ≥2.8                   | without spike, coated |
| PRM-FB-23-180-W-P | 2.3  |  | ≥2.8                   | with spike, coated    |
| PRM-FB-23-200-O-P | 2.3  |  | ≥2.8                   | without spike, coated |
| PRM-FB-23-200-W-P | 2.3  |  | ≥2.8                   | with spike, coated    |
| PRM-FB-23-230-O-P | 2.3  |  | ≥2.8                   | without spike, coated |
| PRM-FB-23-230-W-P | 2.3  |  | ≥2.8                   | with spike, coated    |

## DISPOSABLE BIOPSY FORCEPS ALLIGATOR TEETH

Tissue sampling in digestive and respiratory tracts.



## DISPOSABLE BIOPSY FORCEPS ALLIGATOR TEETH

Alligator teeth design with a larger sampling capacity, no slipping during sampling.

Easier to sample tissues, effectively avoid tearing.


Smooth jaws, effectively decrease the friction to working channel.

### Specifications

| Model                |     | Length (mm) | Endoscope Channel | Characteristics                               |
|----------------------|-----|-------------|-------------------|-----------------------------------------------|
| PRM-FB-18-120-O-O-A1 | 1.8 | 1050        | ≥2.0              | without spike, uncoated, with alligator teeth |
| PRM-FB-18-120-O-P-A1 | 1.8 | 1050        | ≥2.0              | without spike, coated, with alligator teeth   |
| PRM-FB-18-160-O-O-A1 | 1.8 | 1600        | ≥2.0              | without spike, uncoated, with alligator teeth |
| PRM-FB-18-160-O-P-A1 | 1.8 | 1600        | ≥2.0              | without spike, coated, with alligator teeth   |
| PRM-FB-23-160-O-O-A1 | 2.3 | 1600        | ≥2.8              | without spike, uncoated, with alligator teeth |
| PRM-FB-23-160-O-P-A1 | 2.3 | 1600        | ≥2.8              | without spike, coated, with alligator teeth   |
| PRM-FB-23-160-W-O-A1 | 2.3 | 1600        | ≥2.8              | with spike, coated, with alligator teeth      |
| PRM-FB-23-180-O-O-A1 | 2.3 | 1800        | ≥2.8              | without spike, uncoated, with alligator teeth |
| PRM-FB-23-180-O-P-A1 | 2.3 | 1800        | ≥2.8              | without spike, coated, with alligator teeth   |
| PRM-FB-23-180-W-P-A1 | 2.3 | 1800        | ≥2.8              | with spike, coated, with alligator teeth      |
| PRM-FB-23-200-O-O-A1 | 2.3 | 2000        | ≥2.8              | without spike, uncoated, with alligator teeth |
| PRM-FB-23-200-O-P-A1 | 2.3 | 2000        | ≥2.8              | without spike, coated, with alligator teeth   |
| PRM-FB-23-200-W-P-A1 | 2.3 | 2000        | ≥2.8              | with spike, coated, with alligator teeth      |
| PRM-FB-23-230-O-O-A1 | 2.3 | 2300        | ≥2.8              | without spike, uncoated, with alligator teeth |
| PRM-FB-23-230-O-P-A1 | 2.3 | 2300        | ≥2.8              | without spike, coated, with alligator teeth   |
| PRM-FB-23-230-W-P-A1 | 2.3 | 2300        | ≥2.8              | with spike, coated, with alligator teeth      |

## SNARE WITH NET

Grasp and remove foreign bodies from digestive tract.



## SNARE WITH NET



### Specifications

| Model            | O.D.(mm) | Length(mm) | Endoscope Channel(mm) | NET Size (mm) |
|------------------|----------|------------|-----------------------|---------------|
| PRM-FG-18-120-C1 | 1.8      | 1200       | ≥2.0                  | 20*50         |
| PRM-FG-18-120-C2 | 1.8      | 1200       | ≥2.0                  | 25*55         |
| PRM-FG-18-120-C3 | 1.8      | 1200       | ≥2.0                  | 30*60         |
| PRM-FG-18-160-C1 | 1.8      | 1600       | ≥2.0                  | 20*50         |
| PRM-FG-18-160-C2 | 1.8      | 1600       | ≥2.0                  | 25*55         |
| PRM-FG-18-160-C3 | 1.8      | 1600       | ≥2.0                  | 30*60         |
| PRM-FG-25-230-C1 | 2.5      | 2300       | ≥2.8                  | 20*50         |
| PRM-FG-25-230-C2 | 2.5      | 2300       | ≥2.8                  | 25*55         |
| PRM-FG-25-230-C3 | 2.5      | 2300       | ≥2.8                  | 30*60         |
| PRM-FG-25-230-C4 | 2.5      | 2300       | ≥2.8                  | 40*75         |
| PRM-FG-25-230-C5 | 2.5      | 2300       | ≥2.8                  | 40*65         |





## GRASPING FORCEPS

Suit for grasping foreign body in complex shape(e.g.retrieval of non-regular shaped stent.)



## GRASPING FORCEPS



## GRASPING FORCEPS

### Specifications

| Model            | O.D. | Length (mm) | Endoscope Channel (mm) | Characteristics                         |
|------------------|------|-------------|------------------------|-----------------------------------------|
| PRM-FG-18-70-A4  | 1.8  | 70          | ≥2.0                   | A4: Pelican Teeth                       |
| PRM-FG-23-160-A1 | 2.3  | 1600        | ≥2.8                   | A1: Alligator Teeth                     |
| PRM-FG-23-180-A1 | 2.3  | 1800        | ≥2.8                   |                                         |
| PRM-FG-23-230-A1 | 2.3  | 2300        | ≥2.8                   |                                         |
| PRM-FG-23-160-A2 | 2.3  | 1600        | ≥2.8                   | A2: Rat Teeth                           |
| PRM-FG-23-180-A2 | 2.3  | 1800        | ≥2.8                   |                                         |
| PRM-FG-23-230-A2 | 2.3  | 2300        | ≥2.8                   |                                         |
| PRM-FG-23-160-A3 | 2.3  | 1600        | ≥2.8                   | A3: Rat Teeth with Alligator Teeth Jaws |
| PRM-FG-23-180-A3 | 2.3  | 1800        | ≥2.8                   |                                         |
| PRM-FG-23-230-A3 | 2.3  | 2300        | ≥2.8                   |                                         |
| PRM-FG-23-160-A4 | 2.3  | 1600        | ≥2.8                   | A4: Pelican Teeth                       |
| PRM-FG-23-180-A4 | 2.3  | 1800        | ≥2.8                   |                                         |
| PRM-FG-23-230-A4 | 2.3  | 2300        | ≥2.8                   |                                         |
| PRM-FG-23-160-B1 | 2.3  | 1600        | ≥2.8                   | B1: 3 Round Prongs                      |
| PRM-FG-23-230-B1 | 2.3  | 2300        | ≥2.8                   |                                         |
| PRM-FG-23-160-B2 | 2.3  | 1600        | ≥2.8                   | B2: 4 Round Prongs                      |
| PRM-FG-23-230-B2 | 2.3  | 2300        | ≥2.8                   |                                         |
| PRM-FG-23-160-B3 | 2.3  | 1600        | ≥2.8                   | B3: 5 Round Prongs                      |
| PRM-FG-23-230-B3 | 2.3  | 2300        | ≥2.8                   |                                         |



Selamet mh.Vardar Cd.No: 5/35 Ecevit Sinem 2 İş Merkezi Osmangazi BURSA - TURKIYE  
Tel : +90 224 252 02 01 Fax : +90 224 252 03 01  
[www.promostmedikal.com](http://www.promostmedikal.com) [info@promostmedikal.com](mailto:info@promostmedikal.com)