



OBSTETRICS & GYNECOLOGY PRODUCTS

www.promostmedikal.com

DISPOSABLE CERVICAL RIPENING BALLOON

Cervical Ripening Balloon Introduce: The Cervical Ripening Balloon is used for physical dilation of cervical canal prior while the cervix is unfavorable for induction. balloons are inflated with sterilized saline. By using the pressure of both balloon, the cervix will synthesize and release

Introduce:

The Cervical Ripening Balloon is used for physical dilation of cervical canal prior while the cervix is unfavorable for induction. balloons are inflated with sterilized saline. By using the pressure of both balloon, the cervix will synthesize and release prostaglandin, and promote the cervix ripening.

Features

- 1.Shorten the process of birth, reduce the pain of parturition
- 2.Medical grade material ensures the highest level of safety
- 3.Obvious effect on promoting the cervix ripening
- 4.Ripens and dilates the cervix without pharmaceuticals
- 5.Silicone body material eliminates the issue of patient sensitively to latex



PRODUCTS CODE	CATHETER DIAM (Fr)	CATHETER LENGH (cm)	BALLOON VOLUME (ml)
PRM-43.01.10000	18	40	80

NEGATIVE PRESSURE DRAINAGE BALL

Introduce: Used for fetch and collect liquid in vivo clinically.

Content: 1. Silicone ball 2. Volume Marker 3. 100cm drainage tubing 4. Clamp control

Features:

1. Easy to control construction
2. Clear housing designed for bleeding observation



PRODUCT CODE	BALLOON VOLUME(ml)
PRM-43.99.10000	200

POSTPARTUM BALLOON WITH RAPID INSTILATION



Postpartum Balloon with Rapid Instillation

Introduce:

Used for Temporary control or reduction of postpartum uterine bleeding when conservative management is warranted. For while the balloon was inflated with sterilized saline, the postpartum balloon applies pressure to the uterine walls that will help to hemostasis.

Features:

1. Easy to operate and position,
2. safe and efficient way to promote hemostasis
3. Can be used for both vaginal and cesarean delivery
4. Silicone body material eliminates the issue of patient sensitively to latex, 2 pcs 60 ml injector, 180 cm extension tube.

PRODUCT CODE	SIZE (fr)	LENGTH(cm)	BALLOON VOLUME (ml)
PRM-43.01.20000	24	55	500
PRM-R-43.01.20004	24	55	500

LEEP LLETZ ELECTRODES



REF/MODEL	Large Radius LEEP/LLETZ Loop Elektrodes	REF/MODEL	Midium Radius LEEP/LLETZ Loop Elektrodes
PRM-AR 200813	Oval LEEP Elektrod 2.0cm W x 0.8 cm D 13 cm Shaft	PRM-AR 100513	Oval LEEP Elektrod 1.0cm W x 0.5 cm D 13 cm Shaft
PRM-AR 201013	Oval LEEP Elektrod 2.0cm W x 1.0 cm D 13 cm Shaft	PRM-AR 101013	Oval LEEP Elektrod 1.0cm W x 1.0 cm D 13 cm Shaft
PRM-AR 201213	Oval LEEP Elektrod 2.0cm W x 1.2 cm D 13 cm Shaft	PRM-AR 101213	Oval LEEP Elektrod 1.0cm W x 1.2 cm D 13 cm Shaft
PRM- AR 201513	Oval LEEP Elektrod 2.0cm W x 1.5 cm D 13 cm Shaft	PRM-AR 101513	Oval LEEP Elektrod 1.0cm W x 1.5 cm D 13 cm Shaft
PRM-AR 202013	Oval LEEP Elektrod 2.0cm W x 2.0 cm D 13 cm Shaft	PRM-AR 150713	Oval LEEP Elektrod 1.5cm W x 0.7 cm D 13 cm Shaft
PRM-AR 251013	Oval LEEP Elektrod 2.0cm W x 1.0 cm D 13 cm Shaft	PRM-AR 151013	Oval LEEP Elektrod 1.5cm W x 1.0 cm D 13 cm Shaft
PRM-AB 513	Ball LEEP Elektrot 5 mm x 13 cm Shaft	PRM-AR 151213	Oval LEEP Elektrod 1.5cm W x 1.2 cm D 13 cm Shaft
PRM-AB 413	Ball LEEP Elektrot 4 mm x 13 cm Shaft	PRM-AR 151513	Oval LEEP Elektrod 1.5cm W x 1.5 cm D 13 cm Shaft

Selamet mh.Vardar Cd.No: 5/35 Ecevit Sinem 2 İş Merkezi Osmangazi BURSA - TURKIYE
Tel : +90 224 252 02 01 Fax : +90 224 252 03 01
www.promostmedikal.com info@promostmedikal.com